


point designs fabrication options

*Due to the variety and nature of materials used, slight variations in color and design may occur. We cannot guarantee an exact match between similar looking silicone and frame colors.

solid silicone colors



Other colors may be requested. Color choice should be expressed to the Fabrication Team.

silicone patterns

Skin tone




Camo

Other colors*: purple, red, winter




Marble



* Other color combinations may be requested. Certain silicone options will have an upcharge (Camo and other special patterns).

frame color options



frame design options



Frame color and/or patterns may include sports logos, special shapes or designs, etc.
Above are examples of what has been requested on previous projects.

point designs fabrication options

suspension options

Please select one from each column if applicable.

closure

Zipper



Silicone Overlap (Flap)



Open Dorsum



Continuous (No Opening)



strap

Frame Strap
(mounted to frame)



Wrist Strap
(tunneled through silicone)



No Strap



style

Dacron Strap
(default)



Watch Band



Magnetic Watch Band



Quick Lacing System
(additional costs)

